Sakamoto 'Kiyoko' Bathing Care Manikin





The lower back smoothly bends anteriorly and posteriorly. The model can be seated in a wheelchair for simulation with a feeling similar to seating a real patient.

The head can be turned laterally to 30 in both directions as well as upward and downward so that the model can take nearly the same postures as the human body.

Since the lower limbs can be freely moved back and forward or abducted, the legs can be spread sufficiently. Practice of assistance in excretion using a bedpan or a urinal can be performed with a true to life feeling.

Various joints including the shoulders, elbows, wrists, hips, knees and ankles move naturally as those of the human body and enable smooth execution of practice of changing clothes and bedding.







- Can be immersed in water.
- Water spontaneously flows into the model so that it can be dipped into water. It drains quickly when lifted from water.
- Tolerates the use of showers, sponge and detergents.

Height: 160cm Weight: 13kg

Material: Vinyl cholride

Part Number IMPM100-2 Description

Sakamoto 'Kiyoko' Bathing Care Manikin



Mentone Educational Centre 10-12 McWilliam Street Springvale Victoria, 3171, AUSTRALIA.

Ph: (03) 9547 6638 Fax: (03) 9553 4562

International

Ph: +61 3 9553 3234 Fax: +61 3 9553 4562