

Sakamoto Musclemate 2 Intramuscular Injection Simulator



Correct or incorrect injection site is indicated with built-in lamps and a buzzer.



Correct

The green lamp is lit when the injection site and the depth are correct.



Incorrect

The red light is flashed, and the buzzer sounds when the injection site is incorrect.



Incorrect Depth

Both the red and green lamps are flashed and the buzzer sounds when the injection site is correct but too deep.

*The light of the LED lamps is visible through the skin.

An intramuscular injection and palpation simulator.

- Realistic rendition of human skin and underlying muscle.
- Indicates the correct or incorrect injection site including the depth with LED lamps.

Preparation

Palpation

Pinching

Deltoid Muscle

Subcutaneous Injection

