

Sakamoto Buttockmate IM Injection Simulator



Practice gluteus injection skill by palpation. (Shibusanbu Method) Upper Outer quadrant of Buttocks

- Facilitation of training in pairs, undertaking the role of nurse and patient alternately.
- Facilitates palpation of iliac crest, greater trochanter and iliac spine.
- Realistic skin and skeletal structure.
- Incorrect injection indicated by red light and buzzer.
- Correct injection technique indicated by green light.
- Unique material means that no needle mark is left in the skin.




Correct injection site by palpation



Fluid can be injected



Use as a stand alone unit



Transparent overlay indicates
correct injection site.



Part Number
IMPM153

Description
Sakamoto Buttockmate IM Injection Simulator

Mentone Educational Centre
10-12 McWilliam Street Springvale
Victoria, 3171, AUSTRALIA.
Ph: (03) 9547 6638 Fax: (03) 9553 4562
International
Ph: +61 3 9553 3234 Fax: +61 3 9553 4562