

Sakamoto Suction Simulator



MENTONE
EDUCATIONAL
ANATOMICAL MODELS AND EDUCATIONAL SUPPLIES




This model is a simulator with a lifelike texture.

Various suction methods can be practised.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

The head can be moved to the left or right



Includes:

- 1 x Suction Simulator
- 1 x Baby Powder
- 1 x Lubricant Oil
- 1 x Container for Making Artificial Sputum
- 1 x Agitator Stick for Mixing Artificial Sputum
- 1 x Syringe with Tube
- 1 x Silicone Cap
- 1 x Flushing Funnel
- 1 x Flushing Tube
- 1 x Tray



MENTONE
EDUCATIONAL
ANATOMICAL MODELS AND EDUCATIONAL SUPPLIES

Height: W68cm x D34cm x H32cm
Weight: 15kg

Part Number	Description
IMPM175	Sakamoto Suction Simulator
IMPM	Artificial Sputum

Mentone Educational Centre
10-12 McWilliam Street Springvale
Victoria, 3171, AUSTRALIA.
Ph: (03) 9547 6638 Fax: (03) 9553 4562
International
Ph: +61 3 9553 3234 Fax: +61 3 9553 4562

Sakamoto Suction Simulator


Various Types of Suction Training are Possible



Suction of secretions in the oral cavity



Suction of secretions in inserted tubes



Suction of secretions in the nasal cavity



Suction of secretions in the tracheal cannula



Washing of the bronchi



Washing of the oral cavity