

Sakamoto Baby



- The limbs have seamless structures.
- The head moves in every direction. (Please hold the model by supporting the head as you would do in holding a neonate.)
- The weight of the head has the same proportion as that of a neonate.
- The anterior fontanelle and ribs can be palpated.
- Water infused into the mouth flows out from the urethra so that changing diapers can be practiced.
- Realistic newborn manikin with weak neck muscles


1. Bathing



2. Baby exercise and feeding



3. Dressing / undressing and changing diapers



4. Anal thermometry and various measurements



Height: 49cm
Weight: 3kg
Material: Silicone

Part Number
IMPM107-3
IMPM107-4

Description
Sakamoto Baby Male Version
Sakamoto Baby Female Version



Mentone Educational Centre
10-12 McWilliam Street Springvale
Victoria, 3171, AUSTRALIA.
Ph: (03) 9547 6638 Fax: (03) 9553 4562
International
Ph: +61 3 9553 3234 Fax: +61 3 9553 4562