Ligature Training Dummy

HOSPITAL, POLICE, PRISON, MENTAL HEALTH







Equip staff with the right training to safely and sensitively manage hanging scenarios when they occur, including the use of a hook knife to properly release a ligature.

Constructed from the same flame retardant material used in Police ballistic/stab vests, offering 4x the strength and abrasion resistance of 16oz canvas or P.V.C. (conforms to OEKO-TEX®) Standard 100, BS EN ISO 12947-2 & 14116: 2008) High density polypropylene webbing (25mm rated to 207 Kg, 50mm rated to 440Kg).

- This manikin has anatomical landmarks including the sternum, rib cage, and substernal notch. The trainer also contains teeth, uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, oesophagus, inflatable lungs and stomach.
- Allows for placement of a nasogastric tube and accepts Combitube®, E.O.A., E.G.T.A., P.T.L. and i-gel® airways.
- Realistic knee joints promoting correct manual handling techniques including preservation of the C-spine.
- Special D-ring feature enables easy suspension without damaging the manikin.



Model	Description	Height	Weight	
RLNLAW	Ligature Dummy	1.8m	25kg	